



Using EncaptureMD Reporting with SHOWCASE®

This guide provides some basic information about running ShowCase integrated with the EncaptureMD clinical reporting system.

Integration Features

- Display a study in ShowCase and the corresponding report displays in EncaptureMD.
- Bring up an EncaptureMD report and the corresponding images display in ShowCase.
- Automatic transfer of ultrasound machine measurement data from ShowCase to EncaptureMD -- for qualified ultrasound machines that transfer measurements in industry standard DICOM SR format.
- Select key images in ShowCase and send them to the EncaptureMD report by using the Edit menu option **Copy Image to Report**.




We recommend a dual monitor configuration so that you can review images on one screen with ShowCase while entering your interpretation into EncaptureMD on the other screen.

EncaptureMD is web based, so requires internet access on the PC.

Configuration

Licensing: If you add EncaptureMD to an existing ShowCase system you will need to re-activate your ShowCase license with the reporting option turned on. (30-day trial is available.)

Source directories: Use the ShowCase Edit Menu **Report Preferences** option to designate directories to use to look for studies matching the EncaptureMD reports. If you have a local ShowCase Image Center, you can set the Image Center as your study source.



If you do not have a ShowCase Image Center, point to the image directories that have your recent studies.




Once you have designated these study sources the integration with EncaptureMD should be automatic with no further configuration required.

Troubleshooting

Problems with communication between the applications

- Make sure that your internet connection is working so that EncaptureMD can run.
- Check the ShowCase Help Menu option “ShowCase Trace Log...” to see if errors have been logged that might help diagnose the problem.
- If communication between the programs does not seem to be working, try restarting ShowCase
- Check your ShowCase license (File Menu, **License Management** option) to make sure your trial has not run out and that a valid license with the reporting option is activated. The top of the license management screen should look like this:



- If none of the above helps, reboot your machine and start both applications again.

Images matched with wrong report

If you find that an image study is bringing up the wrong report or a report bringing up the wrong images,

- Make sure you are entering the correct Patient ID and Study Date in both ShowCase and EncaptureMD. This information is critical for matching the images and reports.
- If somehow a report and image study become paired incorrectly, use the instructions in your EncaptureMD manual for removing the association.

Images not copying to the report

If you are using the Edit Menu, **Copy to Report** option, make sure that you have selected (clicked on) an image. A selected image has a yellow border around it.

SHOWCASE® is a registered trademark of Trillium Technology, Inc.